

**NOTE: This information has been superseded by a newer version of Skyline. For current documentation, run SkylineRunner.exe or SkylineCmd.exe without any command-line arguments. The same information can be found in the Skyline user interface at Help > Documentation > Command Line.**

## Skyline Command-Line Interface

---

The command-line interface for Skyline is called SkylineRunner.exe. It is a tiny shim executable less than 10 KB in size. It requires a full Skyline installation on the computer on which it is run. SkylineRunner simply starts Skyline running without any user interface, pipes the parameter options from the command-line to the running Skyline instance, and prints output from Skyline to the command console. At present, only one instance of SkylineRunner may be executed at a time.

The Skyline instance started by SkylineRunner is independent of any other instances that may already be running on the same machine. It is not necessary to have a visible instance of Skyline running on your computer for SkylineRunner to work.

SkylineRunner is intended for automating tasks, such as quality control, scheduling and refinement, during acquisition. SkylineRunner can open a Skyline document, import a newly acquired data file, and export a report or new method.

### Parameter Options:

The current implementation of SkylineRunner offers the following options:

#### General input/output

--in=path/to/file.sky	Open a Skyline file
--save	Saves any changes to the file
--out=path/to/file.sky	Same as save except writes to the specified file
--batch-commands =path/to/file	Runs a file line by line treating each line like a SkylineRunner input. Useful for automating the execution of multiple commands. The open Skyline file remains active through all commands.
--dir=path/to/folder	Used to specify a default root directory for all other path arguments, other than the directory in which the command is run.
--timestamp	All logging output will be preceded by a time and date.

*Until the section titled Settings Customization all other command line parameters rely on the "in" parameter because they all rely on having a Skyline document open.*

#### Importing results replicates

--import-file=path/to/file	Attach a replicate to the open document
--import-replicate-name=<name>	Name to give the new replicate in an --import-file

	operation.
--import-optimizing=<ce   dp>	Indicates the data being imported contains extra transitions for detecting optimal collision energy or declustering potential.
--import-append	Append the import-file to the given replicate. This is an intention check in case the document already has a replicate with the given name. By default this is set to false. This option only works with the --import-file option.
--import-all=path/to/folder	Imports from a folder all files or sub-folders which are not already in the document, naming each with the base-name of the file or sub-folder, unless the --import-naming-pattern parameter is also supplied.
--import-naming-pattern=reg-ex	A regular expression from which the first group will be used to name replicates in an --import-all operation (e.g. [^_](.*) for everything after the first underscore).
--import-before=<date>	When importing from a folder, only import from files with modified time before the given date.
--import-on-or-after=<date>	When importing from a folder, only import from files with modified time after the given date.
--import-no-join	Import results files to individual .skyd files without joining them to the main document .skyd file. This is useful for distributed processing, as on HPC cluster.
--import-lockmass-positive	Waters lockmass correction $m/z$ for positive ion scans.
--import-lockmass-negative	Waters lockmass correction $m/z$ for negative ion scans.
--import-lockmass-tolerance	Waters lockmass correction tolerance (Da).

### Reintegrate with advanced peak picking models

--reintegrate-model-name=<name>	The name of a scoring model to use for the reintegrate operation. The model can either be pre-defined (e.g. using the Edit > Refine > Reintegrate form) or created automatically during this operation by using --reintegrate-create-model.
--reintegrate-create-model	This option will cause a new model to be created, using the mProphet algorithm with all

	available scores for the results found in the document. (requires --reintegrate-model-name)
--reintegrate-annotate-scoring	Peaks will be annotated with q value and score annotations. (requires --reintegrate-model-name)
--reintegrate-overwrite-peaks	Existing manually integrated peaks will be overwritten with peaks chosen by the reintegration model. (requires --reintegrate-model-name)

### Removing results replicates

--remove-before=<date>	Remove all results from the open document with an acquired time before the given date.
--remove-all	Remove all results from the open document.

### Importing FASTA files

--import-fasta=path/to/file	Import a FASTA file into the open document.
--keep-empty-proteins	Keeps any empty proteins in the open document after importing a FASTA file.

### Importing peptide searches

--import-search-file=path/to/file	Import a peptide search results file into the open document, building a document-specific spectral library. This may be specified multiple times for multiple files. Use --import-fasta argument to add matched peptides as targets.
--import-search-cutoff-score=<cutoff>	Defines a cutoff score (between 0 and 1) to be used when building a spectral library from peptide search results files, where 1 is for highest confidence matches and 0 includes everything. [default 0.95]
--import-search-add-mods	Adds all modifications found in peptide search results files to the open document.

### Adding spectral libraries

--add-library-path=path/to/file	Specify a spectral library to be added to the open document.
--add-library-name=<name>	Name to give the spectral library in an --add-library-path operation.

### Exporting reports

--report-name=<name>	The name of a report to export as it appears in the Skyline Export Report form
--report-file=path/to/file.csv	The path to export the report to. Required if --

	report-name is specified.
--report-format=<CSV   TSV>	CSV for comma-separated reports (or semicolon separated, depending on your localization) or TSV for tab separated reports [default CSV]
--report-invariant	Exports the report with the “Invariant” language setting, using English (US) number formats and header text without spaces, ideal for use with the R statistical programming environment.

### Exporting chromatograms

--chromatogram-file=path/to/file.tsv	The path to the tab delimited file where to export chromatograms.
--chromatogram-precursors	Export precursor ion chromatograms.
--chromatogram-products	Export product ion chromatograms.
--chromatogram-base-peaks	Export base peak chromatograms.
--chromatogram-tics	Export total ion current chromatograms.

### Exporting transition lists

--exp-translist-instrument=<AB Sciex   Agilent   Thermo   Waters>	Export a transition list. This option is required for exporting a transition list and has no default. This option cannot be used with --exp-method-instrument, because you cannot export a method and transition list simultaneously.
---	---

### Vendor-specific transition list options

AB Sciex	--exp-dwell-time=<millis>	Dwell time per transition. This option is required for unscheduled transition lists.
Agilent	--exp-dwell-time=<millis>	Same as above.
Thermo Scientific	--exp-add-energy-ramp	Adds an extra column for energy ramp to the transition list. Optional. Defaults to false.
Waters	--exp-run-length=<minutes>	Run length of the entire gradient in minutes. This option is required for unscheduled experiments.

### Exporting native instrument methods

--exp-method-instrument=<AB SCIEX QTRAP	Export a method. This option is
---	---------------------------------

Agilent 6400 Series   Thermo TSQ   Thermo LTQ   Waters Xevo   Waters Quattro Premier>	required for exporting a method and has no default. This option cannot be used with --exp-translist-instrument, because you cannot export a method and transition list simultaneously.
--exp-template=path/to/file.meth exp dam m	Path of the method template. This can be a file or a directory depending on your instrument. This option is required for method export.

### Vendor-specific method options

AB Sciex Qtrap	--exp-dwell-time=<millis>	Dwell time per transition. This option is required for standard (unscheduled) methods.
Agilent (all instruments)	--exp-dwell-time=<millis>	Same as above
Thermo (all but LTQ)	--exp-run-length=<minutes>	Run length of the entire gradient in minutes. This option is required for unscheduled experiments.
Waters (all instruments)	--exp-run-length=<minutes>	Same as run length above

### Method and transition list options

--exp-file=path/to/file	Path to the method or transition list file (or directory) to export to. This option is required for method and transition list export.
--exp-strategy=<single   protein   buckets>	Strategy for dividing a method into injections. The default is "single".
--exp-method-type=<standard   scheduled   triggered>	Sets a standard, scheduled or triggered method. The default is "standard".
--exp-max-trans=<number>	Maximum number of transitions per injection for export strategies "protein" and "buckets" OR maximum number of simultaneous transitions for scheduled methods. The default is 100.
--exp-optimizing=<ce   dp>	Export a method with extra transitions for finding optimal collision energy or

	declustering potential.
--exp-scheduling-replicate=<name>	Use this only if creating a scheduled or triggered method. The default is to schedule based on an average of all replicates, but if you specify one, the method will be scheduled based on that replicate.
--exp-ignore-proteins	Ignore protein boundaries in creating methods.
--exp-primary-count=<number>	For --exp-method-type=triggered specifies the number of transitions to make primary.

### Publishing to Panorama

--panorama-server=<server url>	URL of the Panorama server to which the --in file is to be published. The URL should contain the protocol (http or https), hostname and port, if required. Examples: <a href="https://panoramaweb.org">https://panoramaweb.org</a> or http://localhost:8080
--panorama-username=<username>	The username/email address for a user with access to publish to the Panorama server.
--panorama-password=<password>	The password for a user with access to publish to the Panorama server.
--panorama-folder=path/to/folder	The path to a folder on the Panorama server to which the file is to be published (e.g. MyProject/MyFolder).

*If the parameters above are used along with parameters to import results files into the document (--import-file or --import-all) the Skyline document will be uploaded to the given Panorama server only if new results are added to the document.*

### Settings Customization

*The below commands do not rely on the "in" parameter because they modify the user settings that are independent of a specific Skyline document.*

--full-scan-precursor-res=<resolving power>	Resolving power of the precursor mass analyzer.
---	---

<code>--full-scan-precursor-res-mz=&lt;m/z value&gt;</code>	The m/z value at which the precursor mass analyzer resolving power is specified. (applies only to orbitrap and ft_icr mass analyzers)
<code>--full-scan-product_res=&lt;resolving power&gt;</code>	Resolving power of the product mass analyzer.
<code>--full-scan-precursor-res-mz=&lt;m/z value&gt;</code>	The m/z value at which the product mass analyzer resolving power is specified. (applies only to orbitrap and ft_icr mass analyzers)
<code>--full-scan-rt-filter-tolerance=&lt;minutes&gt;</code>	The number of minutes on either side of the predicted time or MS/MS IDs, i.e. $\pm$ minutes. Defaults to.
<code>--tool-arguments="&lt;arguments&gt;"</code>	Optional command-line arguments for the tool to be added, used when the tool is executed. (Not applicable to web URL commands)
<code>--tool-initial-dir=path/to/dir</code>	Optional initial directory for the tool to be added, used when the tool is executed. (Not applicable to web URL commands)
<code>--tool-conflict-resolution=&lt;overwrite   skip&gt;</code>	Tells the SkylineRunner how to resolve a tool name conflict, by either overwriting an existing installation or skipping installation of the new tool.
<code>--tool-report=&lt;report-name&gt;</code>	The name of a report in the settings to use as the input report for the tool.
<code>--tool-output-to-immediate-window</code>	When present the tool output is piped to the Immediate Window at runtime.
<code>--report-add=path/to/file.skyr</code>	Adds the report formats from a skyr file. If there are name conflicts the <code>--report-conflict-resolution</code> parameter is required.
<code>--report-conflict-resolution=&lt;overwrite   skip&gt;</code>	Tells the SkylineRunner how to resolve a report name conflict, by either overwriting the existing report or skipping adding the new report.
<code>--tool-add-zip=path/to/file.zip</code>	Import tools from a tool installation ZIP file.
<code>--tool-zip-conflict-resolution=&lt;overwrite  </code>	Specify whether tool conflicts from the provided ZIP file should be resolved by

parallel>	overwriting or installing in parallel. This is for conflicts related to tool versioning and report names.
--tool-zip-overwrite-annotations=<true   false>	Specify whether conflicting custom annotations from the provided ZIP file should overwrite (true) existing annotations or be skipped (false).
--tool-program-macro=<programTitle> Or --tool-program-macro=<programTitle>,<programVersion> Eg. --tool-program-macro=R,2.15.2	Specifies a program title and version to use with the --tool-program-path command. Together these commands are for importing tools from a ZIP file that use the \$(ProgramPath()) macro as their command. For more information see the documentation on External Tools.
--tool-program-path=path/to/file	Specifies the path to an executable on the local machine for the program title and version specified by the --tool-program-macro flag.
--tool-ignore-required-packages	Ignore required packages when installing a tool from a ZIP file.