

Skyline Small Molecule Method Development and CE Optimization

The Skyline Targeted Proteomics Environment provides informative visual displays of the raw mass spectrometer data you import into your Skyline documents. Originally developed for proteomics use, Skyline has been extended to work with generalized molecules. This tutorial explores using Skyline for targeted quantification of small molecules.

In this tutorial, you will learn about development of a multiplexed method for selected energy metabolites on LC-MS/MS (Triple Quad). In the analysis of this dataset you will learn about:

- Targeted Quantification Workflows based on Triple Quad MS (TQMS)
- Starting from a published transition list of putative molecules of interest, collected on a TQMS from a different vendor
- Building an unscheduled Skyline method using collision energy (CE) values from that list
- Using that unscheduled method to create a scheduled method with optimized retention times
- Using that scheduled method to create a final scheduled method with optimized CE values

You may also wish to view the second half of Skyline webinar 16, on which this tutorial is based:

<https://skyline.ms/project/home/software/Skyline/events/2017%20Webinars/Webinar%2016/begin.viaw?>

Skyline aims to provide a vendor-neutral platform for targeted quantitative mass spec research. It can import raw data from the instrument vendors Agilent, SCIEX, Bruker, Shimadzu, Thermo-Scientific and Waters. The ability to import data across various instrument platforms greatly facilitates cross-instrument comparisons and large multi-site studies. This remains equally true in using it to target small molecules, as it has been for years in the field of proteomics.

If you have not already looked at the “Skyline Small Molecule Targets” tutorial you should do so now, in order to pick up a few basics about how Skyline works with small molecule descriptions including chemical formulas and adducts.

Getting Started

To start this tutorial, download the following ZIP file:

<https://skyline.gs.washington.edu/tutorials/SmallMoleculeMethodDevAndCEOptimization.zip>

Extract the files in it to a folder on your computer, like:

C:\Users\bspratt\Documents


This will create a new folder:

C:\Users\bspratt\Documents\SmallMoleculeMethodDevAndCEOptimization


It will contain all the files necessary for this tutorial.

Now:


- Start Skyline. You should see the Start page. It will look something like this:



- Set the user interface control in the upper right corner of the Start Page to “Molecule interface”. This turns off all the proteomics menus and controls that we won’t be using.



Now the start page is not cluttered with proteomics-related controls that we don’t need.



- Click on “Blank Document”.

Importing a Small Molecule Transition List into a Skyline Document

The transition lists in this tutorial are from a high speed HILIC method, based on Guder et al, [Anal Chem](#), 2017 Feb 7;89(3):1624-1631. The raw data for the tutorial itself was then collected on an Acquity UPLC coupled to a Waters Xevo TQ-S triple quad.

This is the published transition list:

Functional Category	Metabolite	KEGG id	Q1 -12C	Q3-12C	Q1-13C	Q3-13C	CE	Mode
Central Metabolism	a-Ketoglutaric acid	C00026	145	101	150	105	5	Negative
Central Metabolism	Phosphoenolpyruvate	C00074	167	79	170	79	29	Negative
Central Metabolism	Dihydroxyacetone-P	C00111	169	79	172	79	29	Negative
Central Metabolism	Pentose-P	C00199	229	79	234	79	45	Negative
Central Metabolism	Hexose-P	C01094	259	79	265	79	53	Negative
Central Metabolism	Seduheptulose 7-P	C05382	289	97	296	97	17	Negative
Central Metabolism	Fructose-1,6-Bisphosphate	C00354	339	241	345	247	16	Negative
Central Metabolism	UDP-N-acetyl-D-Glucosamine	C00043	606	385	623	394	29	Negative
Central Metabolism	Acetyl-CoA	C00024	808	408	831	418	37	Negative
Cofactor metabolism	NAD	C00003	662	540	683	555	20	Negative
Cofactor metabolism	NADP	C00006	742	620	763	635	20	Negative
Nucleotide metabolism	Orotate	C00295	155	111	160	115	9	Negative
Nucleotide metabolism	Dihydroorotate	C00337	157	113	162	117	5	Negative
Nucleotide metabolism	UDP	C00015	403	79	412	79	69	Negative
Amino acid metabolism	GABA	C00334	104	69	108	73	37	Positive
Amino acid metabolism	Phenylpyruvic acid	C00166	165	95	174	101	13	Positive
Amino acid metabolism	Diaminopimelic acid	C00666	191	128	198	134	13	Positive
Central Metabolism	D-Alanyl-Alanine	C00993	161	44	167	46	13	Positive
Cofactor metabolism	D-Pantothenic acid	C00864	220	90	229	93	13	Positive
Cofactor metabolism	Oxidized glutathione	C00127	613	355	633	365	25	Positive
Nucleotide metabolism	Hypoxanthine	C00262	137	55	142	57	37	Positive
Nucleotide metabolism	Guanine	C00242	152	110	157	114	21	Positive
Nucleotide metabolism	UMP	C00105	325	97	334	102	17	Positive
Nucleotide metabolism	cAMP	C00575	330	136	340	141	29	Positive
Nucleotide metabolism	AMP	C00020	348	136	358	141	21	Positive
Nucleotide metabolism	ADP	C00008	428	136	438	141	37	Positive
Nucleotide metabolism	UTP	C00075	485	97	494	102	21	Positive
Nucleotide metabolism	ATP	C00002	508	136	518	141	37	Positive


It provides information about light (12C) and heavy (13C) precursor and fragment m/z and charge for each molecule. For the purposes of this tutorial, we will work only with the Negative mode entries. The collision energy (CE) values are from a Sciex triple quad, we will use these as a starting point even though we will develop our new method on Waters equipment.

With a little effort in Excel or other spreadsheet editor, this can be reformatted so that Skyline can read it (the heavy and light versions of each transition are expected to be on different lines). The result of this reformatting can be found in the provided **Energy_TransitionList.csv** file in the folder you created at the start of this tutorial.

Because Energy_TransitionList.csv is formatted with column headers that Skyline recognizes, you can bypass the usual **Insert Transition List** window and simply copy everything (include the header row) into Skyline's **Targets** window. To do this, perform the following steps:

- Open Energy_TransitionList.csv in any text editor.
- Select All, then Copy. Make sure you have included the header row.
- In Skyline, click in the **Targets** window, then Paste (**Ctrl+V**).

Your Skyline window should now look something like this (note the molecule symbol in the User Interface control in the upper right – Skyline has turned off all proteomics-specific menus and controls):



NOTE: In this instance, the targets are described with m/z and charge values. Skyline can accept higher level descriptions including chemical formulas and heavy isotope labels etc. Having the chemical formula is especially useful when working with full scan data as it allows Skyline to calculate isotopic distributions, but for SRM data such as this, using m/z and charge is perfectly adequate.

Transition settings

Next we have to make sure Skyline's Transition Settings are correctly set for importing the experimental mass spectrometer results. To do this, perform the following steps:

- On the **Settings** menu, click **Transition Settings**.
- Change settings as necessary to match the following:

The image shows a screenshot of the 'Transition Settings' dialog box in Skyline software. The dialog has a title bar with a close button (X) and a tabbed interface with the following tabs: Prediction, Filter, Library, Instrument, and Full-Scan. The 'Prediction' tab is selected. The settings are as follows:

Setting	Value
Precursor mass:	Monoisotopic
Product ion mass:	Monoisotopic
Collision energy:	Waters Xevo
Declustering potential:	None
Optimization library:	None
Compensation voltage:	None

There is also a checkbox labeled 'Use optimization values when present' which is currently unchecked. At the bottom of the dialog are 'OK' and 'Cancel' buttons.

We will perform a collision energy optimization study on a Waters Xevo as part of method development and optimization, since the collision energies are being used from another instrument (Sciex).

Transition Settings

Prediction Filter Library Instrument Full-Scan

Molecules

Precursor adducts:
[M-]

Fragment adducts:
[M-]

Ion types:
f

Precursor m/z exclusion window:
 m/z

Auto-select all matching transitions

OK Cancel

In this experiment we are interested only in negatively charged transitions, this is indicated in the “Precursor adducts” and “Fragment adducts” fields. Skyline’s adduct descriptions support any kind of ionization description (“[M-H]”, “[M+Na]” etc.), including unknown ionization modes (the “[M-]” used here to indicate “negative mode, charge 1, without any known chemical composition to explain it”). The “f” setting in Ion Types on the Filter tab means we are only interested in fragment ion transitions. If you wanted to measure precursor ions as well, it would be “f,p”.

Transition Settings

Prediction Filter Library Instrument Full-Scan

Min m/z: m/z

Max m/z: m/z

Dynamic min product m/z

Method match tolerance m/z: m/z

Firmware transition limit:

Firmware inclusion limit:

Min time: min

Max time: min

OK Cancel

We will allow a generous match tolerance when reading SRM data and assigning chromatograms to our targets.

Now, save the current Skyline document. To do this:

- From the **Files** menu, click on **Save**.
- When prompted, use the filename "EnergyMet_demo.sky".

Exporting unscheduled methods

In general, if you are running Skyline on a computer with the instrument control software for your brand of mass spec installed, Skyline can export native methods using templates that you specify.

If you do not have Waters MassLynx installed, just skip ahead to the next section “Exporting a Transition List”.

We will export two methods, one of 2 minutes and one of 5 minutes. To begin:

- From the **Files** menu, click on **Export**.
- Then click on **Method** to open the Export Method dialog.
- Adjust the settings until the dialog looks like this:

The screenshot shows the 'Export Method' dialog box with the following settings:

- Instrument type: Waters Xevo TQ
- Single method (selected)
- One method per molecule list (unselected)
- Multiple methods (unselected)
- Ignore molecule lists (unchecked)
- Max transitions per sample injection: (empty)
- Methods: 1
- Optimizing: None
- Method type: Standard
- Run duration (min): 2
- Template file: VerifyETemplate.exp

- Click **OK** and save the method as “EnergyMet_2minutes”.
- Repeat these steps, changing Run duration to 5 and saving as “EnergyMet_5minutes”.

Note: the VerifyETemplate.exp file is in the SmallMoleculeMethodDevAndCEOptimization folder you created at the start of this tutorial.

Exporting a transition list

The steps for exporting a transition list are very similar to that for exporting an unscheduled method, except that there is no template file:

- From the **Files** menu, click on **Export**.
- Then click on **Transition List** to open the Export Transition List dialog.
- Adjust the settings until the dialog looks like this:

The screenshot shows the 'Export Method' dialog box with the following settings:

- Instrument type:** Waters Xevo TQ
- Buttons:** OK (highlighted), Cancel
- Method Selection:**
 - Single method
 - One method per molecule list
 - Multiple methods Ignore molecule lists
- Max transitions per sample injection:** (empty text box)
- Methods:** 1
- Optimizing:** None
- Method type:** Standard
- Run duration (min):** 2
- Template file:** VerifyETemplate.exp

- Click **OK** and save the list as “EnergyTL_2minutes.csv”.
- Repeat these steps, changing Run duration to 5 minutes and saving as “EnergyTL_5minutes.csv”.

Importing mass spectrometer runs

At this point, we would collect data for the samples using our 2 and 5 minute gradients. The samples used here are from a commercially available kit and are 1:1 mixes of light:heavy metabolites extracted from E. Coli (Credentialed E. Coli Lysate from Cambridge Isotope Laboratories).


Results of these two sample collections are in a folder named “Unscheduled” in the folder you created at the start of this tutorial. “01a” is a 2 minute gradient, and “02a” is 5 minutes. We will examine them to decide which gradient is most effective for this experiment.

Perform the following steps.

- On the **File** menu, click **Save**. (Ctrl-S)
- On the **File** menu, select **Import** and click on **Results**.

- In the **Import Results** form, choose to import single-injection replicates. For best performance, be sure to select “Many” in the “Files to import simultaneously” control at the bottom of the form. Now click the **OK** button.
- Navigate to the “Unscheduled” folder and select both files.

The **Import Results Files** form should look like:



- Click the **Open** button.


At this point Skyline notices that the replicate names derived from these file names are very similar - they only vary by a single character.

Skyline asks:

The screenshot shows a dialog box titled "Import Results" with a close button (X) in the top right corner. The main text reads: "The files you have chosen have a common prefix and suffix. Would you like to remove some or all of the prefix or suffix to shorten the names used in Skyline?". There are two radio button options: "Do not remove" (unselected) and "Remove" (selected). Below the options are two text input fields: "Common prefix:" containing "ID33140_0" and "Common suffix:" containing "a_WAA253_4814_092017". A text area labeled "Replicate names:" contains the numbers "1" and "2". At the bottom are "OK" and "Cancel" buttons.


- Click **OK** to use the shortened replicate names "1" and "2".

The files should import within 30 seconds or so, leaving your Skyline window looking something like this:



“1” and “2” are not very memorable replicate names, but we can change them to something more descriptive.


- On the **Edit** menu, choose **Manage Results**.
- Click on the first replicate “1”
- Click on the **Rename** button
- Change the name from “1” to “unscheduled_2min”.



Click the **OK** button.

- Repeat these steps for replicate “2”, renaming to “unscheduled_5min”.
- Click the **OK** button in the **Manage Results** dialog.


Now your Skyline window should look like this:



To take advantage of the Skyline summary graphs for viewing individual targets, do the following:

- On the **View** menu, choose **Arrange Graphs** and click **Tiled**.
- On the **View** menu, choose **Peak Areas** and click **Replicate Comparison**.
- On the **View** menu, choose **Retention Times** and click **Replicate Comparison**.
- Click and drag these views over the arrows to dock them above the chromatogram graphs.

The Skyline window should now look something like this:




Comparing Peak Areas and Resolution

Looking at the **Peak Areas – Replicate Comparison** window we can see that the 5 minute gradient has generally greater peak areas than the 2 minute gradient. You may choose to perform technical replicates to be sure this is not random variance. In this case, the 5 minute method seems to show better peak areas for many analytes but not all, likely based on the better separation of the analytes and decreased ion suppression (also called ‘matrix effect’). Based on this, the 5 minute gradient is the probably the better chromatography choice for this experiment.

A closer examination of the chromatographic peaks is also useful:

- Right-click on either of the chromatograms and make sure that **Synchronize Zooming** is enabled.
- Click and drag within either chromatogram to select the time range from 0.8 minutes to 1.8 minutes.

Skyline should look something like this, from which we can see that the 5 minute gradient provides much better peak separation. This is not surprising, and combined with improved signal intensity it justifies the choice of the longer gradient.




Producing a Scheduled Method

Having decided on the 5 minute gradient, we can now produce a scheduled method or transition list. For the purposes of this tutorial we will assume that you do not have the MassLynx instrument control software installed, and will produce a transition list rather than a native method.

First, we want to set the retention time window for export to our transition list.

- From the **Settings** menu, choose **Peptide Settings**
- Select the **Predictions** tab
- Make sure the “**Use measured prediction times when present**” box is checked

- Set the **Time window** value to 1 minute.



The screenshot shows a dialog box titled "Molecule Settings" with a close button (X) in the top right corner. The "Prediction" tab is selected, with other tabs being "Library", "Labels", and "Quantification".

Under the "Prediction" tab, the following settings are visible:


- Retention time predictor:** A dropdown menu is set to "None". To its right is a small calculator icon.
- Use measured retention times when present**
- Time window:** A text input field contains the value "1", followed by the unit "min".
- Ion mobility predictor:** A dropdown menu is set to "None".
- Use spectral library ion mobility values when present**
- Resolving power:** A text input field is currently empty.
- Linear peak width**

At the bottom right of the dialog box, there are two buttons: "OK" and "Cancel".


- Click **OK**

Now we can produce the transition list:

- From the **Files** menu, click on **Export**.
- Then click on **Transition List** to open the Export Transition List dialog.
- Adjust the settings until the dialog looks like this:



- Click **OK**
- You will be asked which replicate to use for timing – we want the 5 minute gradient, so use these settings:



- Click **OK**
- You will be prompted for a file name, save the scheduled transition list as “**scheduled_5min.csv**”.

Collecting Data with the Scheduled Transition List

At this point we would use the newly generated scheduled transition list to collect new results. This has already been done: data was collected for three runs with different light:heavy ratios. There is a 1:1 mix, a 1:2 mix, and a 2:1 mix.

- On the **File** menu, select **Import** and click on **Results**.
- In the **Import Results** form, choose to import single-injection replicates. For best performance, be sure to select “Many” in the “Files to import simultaneously” control at the bottom of the form. Now click the **OK** button.
- Navigate to the “Scheduled” folder and select all three files.
- Click **OK**
- This time when Skyline offers to simplify the replicate names, choose **Do not remove**

The files you have chosen have a common prefix and suffix.
Would you like to remove some or all of the prefix or suffix to shorten the names used in Skyline?

Do not remove
 Remove

Common prefix:


Common suffix:

Replicate names:
ID33140_03a_WAA253_4814_092017
ID33141_03a_WAA253_4814_092017
ID33142_03a_WAA253_4814_092017

- Click **OK**

As we are done with the two minute gradient, we may remove it from the document:


- From the **Edit** menu, choose **Manage Results**
- Select the **unscheduled_2min** replicate and click the **Remove** button



While we are here, some renaming:

- Select the `unscheduled_5min` replicate
- Click the **Rename** button and change the name to "1:1_1" respectively (we're using a naming convention to show the light: heavy mix ratio and run number – though you can use any naming scheme you like)

- Change the others to “1:1_2”, “2:1_2”, and “1:2_2” (note: this renaming can be done in the Document Grid as well).




- Click **OK**

Evaluating the Mass Spec Runs

We can get visual confirmation of the known ratios of the samples, we looking at the **Peak Areas – Replicate Comparison** window.

- In the **Targets** window, click on **Acetyl-CoA**.

We can see that the peak areas have the expected ratios based on the known mix ratios. From the **Retention Times – Replicate Comparison** view we can also see that the chromatography is stable.



Evaluating Linearity

We can evaluate linearity using Skyline's calibration function.

- From the **View** menu, choose **Document Grid**
- In the **Document Grid**, click on the **Reports** control and select **Replicates**
- Set each replicate's **Sample Type** to **Standard**

- Set each replicates **Analyte Concentration** as appropriate. This is expressed as a number rather than a ratio, so set them to **1, 1, 2, 0.5** respectively.


Document Grid: Replicates

Reports ▾ | 1 of 4 | Export... Actions ▾ | Find:

	Replicate	Sample Type	Analyte Concentration
▶	1:1.1	Standard	1
	1:1.2	Standard	1
	2:1.2	Standard	2
	1:2.2	Standard	0.5


We also need to check some quantification settings:

- On the **Settings** menu, click **Peptide Settings**.
- Select the **Quantitation** Tab, and change settings as necessary to match the following:



Now we are ready to examine the calibration curve.

- Close the **Document Grid** view.
- From the **View** menu, select **Calibration Curve**. You should see this:



The three dilution points evaluated here clearly are not as linear as we would like, therefore we would expect there to be some additional work to be done to obtain information on the linear dynamic range of this particular assay.

Collision Energy Optimization

At this point we want to collect collision energy optimization data in hopes of obtaining as much sensitivity as possible for the analytes in the assay; remember these collision energy values were from a completely different instrument vendor. Skyline will help generate a series of scheduled transition lists with collision energies automatically varied around the initial explicit collision energy we set, which we will use to collect new results for our samples. Skyline will then use this data to select the optimal CE values.

We begin by checking our collision energy settings:

- From the **Settings** menu, choose **Transition Settings**
- Select the **Prediction** tab, and click on the **Collision Energy** drop-down control
- Select **Waters Xevo**
- Click on the **Collision Energy** drop-down control again and select **<Edit Current...>**
- This opens the **Edit Collision Energy Equation** window

- Set the **Step Size** to **2** and the **Step Count** to **5**. This will generate a series of methods (or transition list) where the collision energy is iterated up and down in **2** volt increments, for a total of **5** steps in each direction. A generic recommendation is to start with relatively large Step Size (2 or 3 volts) and then potentially repeat the CE optimization with a small Step Size (1 V) in order to perform fine optimization.

Edit Collision Energy Equation [X]

Name:

Regression parameters:

	Charge	Slope	Intercept
▶	2	0.037	-1.066
	3	0.036	-1.328
*			

Optimization: _____

Step size: Step count:

- Click **OK**

- Back in the **Prediction** tab, make sure the **Use optimization values when present** box is checked
- We aren't measuring precursors, so set the **Optimize by** control to **Transitions**

Transition Settings

Prediction Filter Library Instrument Full-Scan

Precursor mass: Monoisotopic

Product ion mass: Monoisotopic

Collision energy: Waters Xevo

Declustering potential: None

Optimization library: None

Compensation voltage: None

Use optimization values when present

Optimize by: Transition

OK Cancel

- Click **OK**

Now we can proceed to exporting the scheduled transition lists. Again, if you are working on a machine with the instrument control software installed, it is generally preferable to export native methods, but for the purposes of this tutorial we will export transition lists.

- From the **File** menu, click on **Export** and choose **Transition List**

- Change settings as needed to match this:

The screenshot shows the 'Export Transition List' dialog box with the following settings:

- Instrument type: Waters
- Single method:
- One method per molecule list:
- Multiple methods: Ignore molecule lists:
- Max concurrent transitions: 100
- Methods: 5
- Optimizing: Collision Energy
- Method type: Scheduled

Note that we are exporting multiple methods in order to limit the number of concurrent transitions. We do this to ensure enough points across each chromatographic peak in Skyline.

- Click **OK**
- You should see the **Scheduling Data** window, set it to use the retention times from the **1:2_2** replicate for scheduling:

The screenshot shows the 'Scheduling Data' dialog box with the following settings:

- Use retention time average:
- Use values from a single data set:
- Replicate: 1:2_2

- Click **OK**
- You will be asked for a file name. Since we are producing multiple transition lists, we actually just want to provide a partial name as a basis for naming the files produced. Set the name as **TL_CE_Opt** and click **Save**
- This should produce five files: TL_CE_Opt_0001.csv, TL_CE_Opt_0002.csv, etc.

Now we take these new transition lists (five in total) and run our 1:1 sample once with each of these methods, and import those results. These results are available in the “CE Optimization” folder of the folder you created at the start of this tutorial:

- From the **File** menu, click **Import** and select **Results**
- We want to combine the five CE optimization runs into a single chromatogram. Adjust the settings to look like this:

Import Results

Add single-injection replicates in files

Add multi-injection replicates in directories

Add one new replicate

Name:
CE Optimization

Optimizing:
Collision Energy

Add files to an existing replicate

Name:

Files to import simultaneously:
Many

Show chromatograms during import


Retry after import failure

OK

Cancel


- Click **OK**

- Select all five files in the **CE Optimization** folder in your tutorial download folder:



- Click **OK**


The files load as a single replicate, and now Skyline should show a new replicate “CE Optimization”.



Let’s investigate the results of the CE optimization:

- In the **Peak Areas – Replicate Comparison** window, **right-click** then choose **Transitions** and enable **Single**
- In that same window, **right-click** again, choose **Transitions**, and enable **Split Graph**


Skyline should look like this:



The **Peak Areas – Replicate Comparison** window now shows us the heavy and light transitions in separate displays. Within the CE Optimization replicate, the individual bars represent the effects of stepping the collision energy. Red represents the original CE value (as used with the Sciex instrument in the literature), the other bars show the effect of 2eV steps away from that. We can see that for Acetyl-CoA, at least, the original value or the -2eV step both provide good peak areas. This is easier to see if we free up room in the graph by switching off the legend:

- In the **Peak Areas – Replicate Comparison** window, **right-click** then disable **Legend**


Skyline should look like this:



Now we can explore the other molecules to see how well we have optimized collision energies. To do this, click on any molecule in the **Targets** window, or use the up and down arrows after the first click to scroll through the list. Let's look at **Pentose-P** in particular:

- In the **Targets** window, click on **Pentose-P**


Looking at the chromatograms for the various CE step values for Pentose-P it's clear that the CE value from the literature isn't optimal for this molecule on this equipment. As the best CE value (**Step-5**, the blue bar) is at the limit of the range we explored, and given the trend toward that end of the range, further investigation may be warranted to arrive at a truly optimal CE for that compound. For that reason, starting with a wider set of steps (such as 3V) may also be warranted when translating between instrument vendors.



Even so, the **Step-5** CE value (blue bar) is a clear improvement over the initial (red bar) value, and we should proceed to creating a new scheduled transition list that uses the most effective observed CE values.

- From the **File** menu, click on **Export** and choose **Transition List**

This time we want just a single method, still scheduled but not optimizing anything this time:




The 'Export Transition List' dialog box is shown with the following settings:

- Instrument type: Waters
- Method selection: Single method, One method per protein, Multiple methods (with Ignore proteins)
- Max concurrent transitions: 100
- Methods: 1
- Optimizing: None
- Method type: Scheduled

Buttons for OK and Cancel are visible.

- Click **OK**
- You should see the **Scheduling Data** window. Use the retention times from the **CE Optimization** replicate for scheduling. The optimal CE values will also be taken from this replicate.



The 'Scheduling Data' dialog box is shown with the following settings:

- Selection: Use retention time average, Use values from a single data set
- Replicate: CE Optimization

Buttons for OK and Cancel are visible.

- Click **OK**
- You will be asked for a file name. Set the name as **TL_CE_Final.csv** and click **Save**

Comparing the Original CE Values with the Optimized CE Values

At this point, it is interesting to compare where we began – using SRM transitions and CE values from the literature for a Sciex mass spec published as a simple flat file; to where we have arrived – a retention time scheduled and CE-optimized SRM method for 18 energy metabolites and their internal standards for the Waters Xevo TQ-S.

- Open the original (Sciex instrument sourced) transition list **Energy_TransitionList.csv** and the final (Waters TQMS optimized) transition list **TL_CE_Final.csv** in Excel or any other suitable viewer
- Compare the CE values for various values. A couple of points of interest:
 - Interestingly, Acetyl-CoA has the same optimal CE for both the Sciex and Waters instruments
 - Pentose-P has quite different values: 45eV for Sciex and 35eV for Waters.

If we did wish to further optimize CE for Pentose-P (recall that the CE value we identified as optimal was at the edge of the tested range), we could repeat the process using the **TL_CE_Final.csv** scheduled transition list as the starting point for a new round of CE optimization.

Notes: As mentioned previously, in general when performing iterative optimization, it's good to start with large step values to test a broad range of CE values, then move to smaller steps in later iterations to narrow in on the final values. When performing CE optimization, if the “**Use Optimization Values**” button is checked under **Settings** menu, **Transition Settings, Prediction**, you can expect that the new exported method or transition list will automatically incorporate the optimum collision energy into the method. No manual curation of the CE optimization data is required, other than examination to see if a wider CE optimization range needs to be explored.

Conclusion

In this tutorial, you have learned how to create a Skyline document that targets stable isotope labeled small molecules specified as only precursor m/z, product ion m/z, and collision energy values, from a literature citation. You performed retention time scheduling and collision energy optimization for small molecules by importing a multi-replicate data set from a Waters Xevo TQ-S using initial CE values from a Sciex triple quad. You learned how many existing Skyline features created initially for targeted proteomics use can now be applied to generalized molecule data. Non-proteomic molecule support is still a relatively new feature area for Skyline. As such, you can expect it to continue improving rapidly.