

Skyline Small Molecule Quantification

The Skyline Targeted Proteomics Environment provides informative visual displays of the raw mass spectrometer data you import into your Skyline documents. Originally developed for proteomics use, Skyline has been extended to work with non-proteomic molecules. This tutorial explores a relatively straightforward example of using Skyline for targeted quantification of a single small molecule using an external calibration curve and stable-isotope labeled internal standard.

In this tutorial, you will learn about Targeted Quantification based on TQ-MS (in this example, out of crashed plasma) starting from a method you may already be running (e.g. a pharmacokinetic assay). In the analysis of this dataset you will learn:

- Insertion of simple set of known transitions
- Data Analysis and peak integration for non-proteomic molecules
- Small Molecule Quantification workflow in Skyline

You may also wish to view the first half of Skyline webinar 16, on which this tutorial is based:

[https://skyline.ms/project/home/software/Skyline/events/2017%20Webinars/Webinar%2016/begin.vi
w?](https://skyline.ms/project/home/software/Skyline/events/2017%20Webinars/Webinar%2016/begin.vi
w?)

Skyline aims to provide a vendor-neutral platform for targeted quantitative mass spec research. It can import raw data from the instrument vendors Agilent, SCIEX, Bruker, Shimadzu, Thermo-Scientific and Waters. The ability to import data across various instrument platforms greatly facilitates cross-instrument comparisons and large multi-site studies. This remains equally true in using it to target small molecules, as it has been for years in the field of proteomics.

If you have not already looked at the “Skyline Small Molecule Targets” tutorial you should do so now, in order to pick up a few basics about how Skyline works with small molecule descriptions including chemical formulas and adducts.

Getting Started

To start this tutorial, download the following ZIP file:

<https://skyline.gs.washington.edu/tutorials/SmallMoleculeQuantification.zip>

Extract the files in it to a folder on your computer, like:

C:\Users\bspratt\Documents


This will create a new folder:

C:\Users\bspratt\Documents\SmallMoleculeQuantification


It will contain all the files necessary for this tutorial.

Now:


- Start Skyline. You should see the Start page. It will look something like this:



- Set the user interface control in the upper right corner of the Start Page to “Molecule interface”.




Now the start page is not cluttered with proteomics-related controls that we don't need.



- Click on “Blank Document”.

Experimental Layout

This experiment was designed according to the FDA Guidance on Bioanalytical Method Validation, and as such contains much more than just the study samples. A full description of the plate layout and run order typically utilized for such a study has been published (<https://www.ncbi.nlm.nih.gov/pubmed/29039849>). Briefly, the samples for this dataset were laid out in a 96 well plate as follows:



Blanks, or 'zero' standards, contain only the internal standard, double blanks contain no standard at all.

Calibration curve samples are a dilution series for calibration.

QC samples are "known unknowns". These are quality control samples which we will treat as unknowns in our study. In actuality, we know what the results should be and thus they can be used to check the accuracy of our measurements.

Serum SPQC is a serum pooled QC, a pooling of all study samples which is run at several points at the start, middle and end of the experiment to verify that quantitative reproducibility is constant across the study.

The NIST SRM 1950 is a pooled plasma standard from the National Institute for Standards and Technology, which is available to all researchers as a reference standard for 'normal' plasma metabolite measurements. It serves as a reference between studies in different laboratories.

Injections were performed in this order:

- Serum Sample Injection Sequence
- double blank
 - PBS 'zero' samples
 - Calibration curve Low to High
 - Low, Mid, High QC Samples
 - Global Reference QC
 - ① 4765 Serum SPQC
 - ① NIST SRM 1950
 - Study Samples 1 to 38
 - Low, Mid, High QC Samples
 - Global Reference QC
 - ② 4765 Serum SPQC
 - Study Samples 39 to 76
 - Global Reference QC
 - ③ 4765 Serum SPQC
 - ③ NIST SRM 1950
 - High, Mid, Low QC Samples
 - Calibration Curve Low to High

In all, 113 injections were used in collecting the mass spec data for these samples.

Internal Standards

For this study we have just two targets: a molecule and an internal standard, which is an isotopically labeled variant of the molecule and thus co-elutes. It is also possible to establish a relationship between unrelated molecules by declaring one of the molecules as a surrogate standard. The surrogate standard method is covered in the "Skyline High Resolution Metabolomics" tutorial.

Importing a Small Molecule Transition List into a Skyline Document

The easiest way to get a small molecule transition list into a Skyline document is to start with an empty document and use the **Edit > Insert > Transition List** menu item.

Transition list insert

To begin, do the following:

- On the Skyline **Edit** menu, choose **Insert** and click **Transition List**.

Skyline will show the **Insert** form. You may have a different column selection and order from previous uses of Skyline:

The screenshot shows the 'Insert' dialog box in Skyline. The title bar reads 'Insert' with a close button. Below the title bar is a tab labeled 'Transition List'. The main area contains a table with the following columns: Molecule List Name, Precursor Name, Precursor Formula, Precursor Adduct, Precursor m/z, Precursor Charge, Product Name, Product Formula, Product Adduct, Product m/z, Product Charge, Label Type, Explicit Retention Time, Explicit Retention Time Window, Explicit Collision Energy, Note, and InChiKey. The first row of the table contains a blue cell with '>>' in the 'Molecule List Name' column, and all other cells are empty. At the bottom of the dialog, there are buttons for 'Columns...', 'Help', 'Check for Errors', 'Insert', and 'Cancel'.

Normally you would copy and paste a transition list from Excel or some other external source, but in this case we have a small enough transition list that we can enter it by hand.

You can see that there are some extra column headers in the **Insert** form, and the column order is not the desired order. Both issues are easy to correct:

- Click the **Columns** button and uncheck the columns that do not appear in the spreadsheet.

This should result in a column picking menu like the one shown below:

<input checked="" type="checkbox"/>	Molecule List Name
<input checked="" type="checkbox"/>	Precursor Name
<input type="checkbox"/>	Precursor Formula
<input type="checkbox"/>	Precursor Adduct
<input checked="" type="checkbox"/>	Precursor m/z
<input checked="" type="checkbox"/>	Precursor Charge
<input type="checkbox"/>	Product Name
<input type="checkbox"/>	Product Formula
<input type="checkbox"/>	Product Adduct
<input checked="" type="checkbox"/>	Product m/z
<input checked="" type="checkbox"/>	Product Charge
<input checked="" type="checkbox"/>	Label Type
<input checked="" type="checkbox"/>	Explicit Retention Time
<input type="checkbox"/>	Explicit Retention Time Window
<input checked="" type="checkbox"/>	Explicit Collision Energy
<input type="checkbox"/>	Note
<input type="checkbox"/>	InChiKey
<input type="checkbox"/>	CAS
<input type="checkbox"/>	HMDB
<input type="checkbox"/>	InChi
<input type="checkbox"/>	SMILES
<input type="checkbox"/>	S-Lens
<input checked="" type="checkbox"/>	Cone Voltage
<input type="checkbox"/>	Explicit Drift Time (msec)
<input type="checkbox"/>	Explicit Drift Time High Energy Offset
<input type="checkbox"/>	Explicit Ion Mobility
<input type="checkbox"/>	Explicit Ion Mobility Units
<input type="checkbox"/>	Explicit Ion Mobility High Energy Offset
<input type="checkbox"/>	Collisional Cross Section (sq A)
<input type="checkbox"/>	Explicit Compensation Voltage
<input type="checkbox"/>	Explicit Declustering Potential

Name

Columns... Help

Next do the following to reorder the columns in the **Insert** form:

- Click and drag each column header you want to move to the order matching the spreadsheet. Once you have selected and arranged your columns, the insert form should now appear as shown below:

Enter the following values in the **Insert** form (or better yet, copy and paste from this PDF):

DrugX, Drug, light, 283.04, 1, 129.96, 1, 26, 16, 2.7

DrugX, Drug, heavy, 286.04, 1, 133.00, 1, 26, 16, 2.7

If you accidentally got the column order wrong, then you will see an error at this point. Otherwise, the **Insert** form should look like this:


- Click the **Insert** button.

NOTE: In this instance, the targets are described with m/z and charge values. Skyline can accept higher level descriptions including chemical formulas and heavy isotope labels etc. Having the chemical formula is especially useful when working with full scan data as it allows Skyline to calculate isotopic distributions, but for SRM data such as this m/z and charge are perfectly adequate.

To see the newly imported targets in full detail:

- From the **Edit** menu, click on **Expand All**
- Then click on **Precursors**.


Your Skyline window should now look like this:



Transition settings

Next we have to make sure Skyline's Transition Settings are correctly set for importing the experimental mass spectrometer results. To do this, perform the following steps:

- On the **Settings** menu, click **Transition Settings**.
- Change settings as necessary to match the following:



Note: the “f” setting in Ion Types on the Filter tab means we are only interested in fragment ion transitions. If you wanted to measure precursor ions as well, it would be “f,p”.

In the Instrument tab, make sure that the minimum and maximum m/z values make sense for your actual instrument. The purpose of these settings is to alert you when you propose to monitor transitions that your mass spectrometer is not actually able to measure.

Method match tolerance is another important setting in the Instrument tab – it determines how closely a measured m/z value must match the transition value to be considered a match.


Now we are ready to import the experimental mass spectrometer results.

Importing mass spectrometer runs

This experiment has 113 mass spec result files associated with it. In cases like this, it can be useful to initially import just a handful of the unknowns along with all of the Calibration Curve and Quality Control (QC) files. This lets us start out with a less complicated document while we verify chromatography and calibration curves etc.

Perform the following steps.

- On the **File** menu, click **Save**. (Ctrl-S)
- Save this document as “SMQuant_v1.sky” in the tutorial folder you created.
- On the **File** menu, select **Import** and click on **Results**.
- In the **Import Results** form, choose to import single-injection replicates. For best performance, be sure to select “Many” in the “Files to import simultaneously” control at the bottom of the form.



Import Results

Add single-injection replicates in files

Optimizing:
None

Add multi-injection replicates in directories

Add one new replicate

Name:
[Text Box]

Add files to an existing replicate

Name:
[Dropdown]

Files to import simultaneously:
Many

Show chromatograms during import


Retry after import failure

OK

Cancel


- Now click the **OK** button.
- You should see a file selection dialog with all the raw data folders in the tutorial folder. Select the last 16 unknowns samples and all the QC samples by clicking the “80_0_1_1_00_1021523383.raw” file and then holding down the Shift key and clicking the last file in the list.

The **Import Results Files** form should look like:



- Click the **Open** button.


The files should import within 30 seconds or so, leaving your Skyline window looking something like this:



To take advantage of the Skyline summary graphs for viewing individual targets, do the following:

- On the **View** menu, choose **Peak Areas** and click **Replicate Comparison**.
- On the **View** menu, choose **Retention Times** and click **Replicate Comparison**.
- Click and drag these views to the arrows to dock them above the chromatogram graphs.
- Select the first target "Drug" in the **Targets** view.


The Skyline window should now look something like this:



Checking Peak Integration

Looking at the **Retention Times – Replicate Comparison** window we can see by the outliers that in a few replicates Skyline has not chosen peaks at the retention times we might have expected.

Clicking on the first outlier in the **Retention Times – Replicate Comparison** window changes the selected replicate to DoubleBlank1, where we would not actually expect Skyline to find a good peak for either the light or heavy form of the drug.




Clicking on the other two outliers reveals that DoubleBlank2 and DoubleBlank3 are the other replicates with apparently poor chromatography. But, of course, as these are double blanks we don't expect any actual peaks in these replicates, so we will manually adjust the integration for each of the double blank replicates to center on the low signal area at 2.7 minutes.

Adjusting Peak Integration

To adjust peak integration, select the DoubleBlank1 replicate then position the mouse below the Retention Time axis (the cursor will change its shape to this: **+**). Click below the Retention Time axis at about 2.65 minutes and drag to about 2.75 minutes. The peak boundaries will change to these new values, and the original boundaries will be marked by a shaded area.


Repeat for the other two DoubleBlank replicates. This should leave Skyline showing something like:



Preparing for Quantitation

Now we can proceed to quantitation. To do this, perform the following steps:

- On the **Settings** menu, click **Molecule Settings**.
- Select the **Quantitation** Tab, and change settings as necessary to match the following:



The screenshot shows the 'Molecule Settings' dialog box with the 'Quantification' tab selected. The settings are as follows:

- Regression fit: Linear
- Normalization method: Ratio to Heavy
- Regression weighting: $1 / (x * x)$
- MS level: All
- Units: uM
- Figures of merit:
 - Max LOQ bias: [] %
 - Max LOQ CV: [] %
 - Calculate LOD by: None

Buttons for 'OK' and 'Cancel' are visible at the bottom right.

We will use a linear regression fit, normalizing to the heavy labeled drug. Skyline offers options for weighting across the curve as a function of x : None, $1/x$, and $1/(x*x)$. We will use a Regression Weighting of $1/(x*x)$. The Units setting is for display purposes, and can be set to any value that makes sense for your experiment. The data experiment was calibrated in micromolar, so we set this to **uM**.

- Click **OK**.

We are still not quite ready to view a calibration curve, however, as we have not yet declared the sample types of the various replicates for Skyline.

Declaring Sample Types for Calibration Curve Display

We will use the Document Grid to examine and add information about the various replicates. The Document Grid is a highly useful tool in Skyline, providing spreadsheet-like views of many document details, much of which can be edited right there in the grid. In this case, we are interested in details of the various replicates, so:

- On the **View** menu, select **Document Grid**
- Click on the **Reports** control in the upper left corner of the grid, and select **Replicates**

The document grid should look like this:

	Replicate	Sample Type	Analyte Concentration
▶	80_0_1_1_00_1...	Unknown	
	81_0_1_1_00_1...	Unknown	
	82_0_1_1_00_1...	Unknown	
	83_0_1_1_00_1...	Unknown	
	84_0_1_1_00_1...	Unknown	
	87_0_1_1_00_1...	Unknown	
	88_0_1_1_00_1...	Unknown	
	89_0_1_1_00_1...	Unknown	
	90_0_1_1_00_1...	Unknown	
	91_0_1_1_00_1...	Unknown	
	92_0_1_1_00_1...	Unknown	
	93_0_1_1_00_1...	Unknown	
	94_0_1_1_00_1...	Unknown	
	95_0_1_1_00_1...	Unknown	
	96_0_1_1_00_1...	Unknown	
	Blank_01	Unknown	

- Expand the Document Grid if needed so that you can see all the replicates at one.
- If needed, alphabetize the list by clicking on the Replicate column header and choose “Sort Ascending”.

Now we can set the various sample types. The replicates with the unknowns are ready to go by default.

- Click on the Sample Type column for “Blank_01”. Change the sample type from “Unknown” to “Blank”.
- Now shift+click on the Sample Type column for “Blank_03” to select all three Blank replicates at once.
- Right-click and select **Fill Down**. Everything in the multiple selection now has the same value as the first item in the selection.

Repeat as needed (or, skip ahead to the table below):

- Set the “Cal_” replicates as **Standard**
- Set the “DoubleBlank_” replicates as **Double Blank**
- Set the “QC_” replicates as **Quality Control**

Recall that the SPQC_ replicates are quality control in a different sense (a pooling of all study samples), so we leave those as **Unknown**.

Analyte concentrations can be entered by hand, but it is much easier to copy and paste them into the grid. Locate the **Concentrations.xlsx** file in the directory you created for this tutorial and open it in Excel. It should look something like this:

Blank_01	Blank	
Blank_02	Blank	
Blank_03	Blank	
Cal_1_01	Standard	10
Cal_1_02	Standard	10
Cal_2_01	Standard	20
Cal_2_02	Standard	20
Cal_3_01	Standard	100
Cal_3_02	Standard	100
Cal_4_01	Standard	200
Cal_4_02	Standard	200
Cal_5_01	Standard	400
Cal_5_02	Standard	400
Cal_6_01	Standard	600
Cal_6_02	Standard	600
Cal_7_01	Standard	800
Cal_7_02	Standard	800
DoubleBlank1	Double Blank	
DoubleBlank2	Double Blank	
DoubleBlank3	Double Blank	
QC_High_01	Quality Control	589
QC_High_02	Quality Control	589
QC_High_03	Quality Control	589
QC_Low_01	Quality Control	121
QC_Low_02	Quality Control	121
QC_Low_03	Quality Control	121
QC_Mid_01	Quality Control	346
QC_Mid_02	Quality Control	346
QC_Mid_03	Quality Control	346
SPQC_01	Unknown	
SPQC_02	Unknown	
SPQC_03	Unknown	

- Select the contents of the table (cells A1 to C32)
- Make sure the order matches the Document Grid
- Click on the “Blank_01” cell and paste

The document grid should look like this when you are done:

Document Grid: Replicates


Reports ▾ | | | 47 of 47 | Export... | Actions ▾ | Find:

Replicate	Sample Type	Analyte Concentration
80 0 1 1 00 1...	Unknown	
81 0 1 1 00 1...	Unknown	
82 0 1 1 00 1...	Unknown	
83 0 1 1 00 1...	Unknown	
84 0 1 1 00 1...	Unknown	
87 0 1 1 00 1...	Unknown	
88 0 1 1 00 1...	Unknown	
89 0 1 1 00 1...	Unknown	
90 0 1 1 00 1...	Unknown	
91 0 1 1 00 1...	Unknown	
92 0 1 1 00 1...	Unknown	
93 0 1 1 00 1...	Unknown	
94 0 1 1 00 1...	Unknown	
95 0 1 1 00 1...	Unknown	
96 0 1 1 00 1...	Unknown	
Blank_01	Blank	
Blank_02	Blank	
Blank_03	Blank	
Cal_1_01	Standard	10
Cal_1_02	Standard	10
Cal_2_01	Standard	20
Cal_2_02	Standard	20
Cal_3_01	Standard	100
Cal_3_02	Standard	100
Cal_4_01	Standard	200
Cal_4_02	Standard	200
Cal_5_01	Standard	400
Cal_5_02	Standard	400
Cal_6_01	Standard	600
Cal_6_02	Standard	600
Cal_7_01	Standard	800
Cal_7_02	Standard	800
DoubleBlank1	Double Blank	
DoubleBlank2	Double Blank	
DoubleBlank3	Double Blank	
QC_High_01	Quality Control	589
QC_High_02	Quality Control	589
QC_High_03	Quality Control	589
QC_Low_01	Quality Control	121
QC_Low_02	Quality Control	121
QC_Low_03	Quality Control	121
QC_Mid_01	Quality Control	346
QC_Mid_02	Quality Control	346
QC_Mid_03	Quality Control	346
SPQC_01	Unknown	

Inspecting the Calibration Curve

Now we are ready to examine the calibration curve.


- Close the **Document Grid** view.
- From the **View** menu, select **Calibration Curve**. You should see this:



The note about the selected replicate missing transitions is expected, as the currently selected replicate is a double blank.

Looking at the graph, we note that the Unknowns loaded so far have relatively low dynamic range.

We also note that some of the calibration samples look like they might be outliers. We can use the Document Grid to get a qualitative sense of that, and to exclude any samples that are not suitable:

- On the **View** menu, select **Document Grid**
- Click on the **Reports** control in the upper left corner of the grid, and select **Replicates**
- Click again on the **Reports** control in the upper left corner of the grid, and select **Customize Report**
- Click on the search button  and type **accuracy** into the **Find What** field
- Click **Find Next**
- In the **Customize Report** form, **Accuracy** should be highlighted under Quantification
- Click on the checkbox to enable **Accuracy**

- Just above Quantification, click on the checkbox to enable **Exclude from Calibration**
- In the **View Name** field at the top of the **Customize View** form, type **Replicates_custom_quant**
- Click **OK**

The Document Grid view should now look like this:

Document Grid: Replicates_custom_quant


Reports | 47 of 47 | Export... | Actions | Find:

Replicate	Sample Type	Analyte Concentration	Accuracy	Exclude From Calibration
80 0 1 1 00 1...	Unknown		#N/A	<input type="checkbox"/>
81 0 1 1 00 1...	Unknown		#N/A	<input type="checkbox"/>
82 0 1 1 00 1...	Unknown		#N/A	<input type="checkbox"/>
83 0 1 1 00 1...	Unknown		#N/A	<input type="checkbox"/>
84 0 1 1 00 1...	Unknown		#N/A	<input type="checkbox"/>
87 0 1 1 00 1...	Unknown		#N/A	<input type="checkbox"/>
88 0 1 1 00 1...	Unknown		#N/A	<input type="checkbox"/>
89 0 1 1 00 1...	Unknown		#N/A	<input type="checkbox"/>
90 0 1 1 00 1...	Unknown		#N/A	<input type="checkbox"/>
91 0 1 1 00 1...	Unknown		#N/A	<input type="checkbox"/>
92 0 1 1 00 1...	Unknown		#N/A	<input type="checkbox"/>
93 0 1 1 00 1...	Unknown		#N/A	<input type="checkbox"/>
94 0 1 1 00 1...	Unknown		#N/A	<input type="checkbox"/>
95 0 1 1 00 1...	Unknown		#N/A	<input type="checkbox"/>
96 0 1 1 00 1...	Unknown		#N/A	<input type="checkbox"/>
Blank_01	Blank		#N/A	<input type="checkbox"/>
Blank_02	Blank		#N/A	<input type="checkbox"/>
Blank_03	Blank		#N/A	<input type="checkbox"/>
Cal_1_01	Standard	10	103%	<input type="checkbox"/>
Cal_1_02	Standard	10	97.4%	<input type="checkbox"/>
Cal_2_01	Standard	20	95.7%	<input type="checkbox"/>
Cal_2_02	Standard	20	100.5%	<input type="checkbox"/>
Cal_3_01	Standard	100	114.6%	<input type="checkbox"/>
Cal_3_02	Standard	100	106.4%	<input type="checkbox"/>
Cal_4_01	Standard	200	105.1%	<input type="checkbox"/>
Cal_4_02	Standard	200	101.8%	<input type="checkbox"/>
Cal_5_01	Standard	400	73.4%	<input type="checkbox"/>
Cal_5_02	Standard	400	76.6%	<input type="checkbox"/>
Cal_6_01	Standard	600	108.9%	<input type="checkbox"/>
Cal_6_02	Standard	600	112.8%	<input type="checkbox"/>
Cal_7_01	Standard	800	107.5%	<input type="checkbox"/>
Cal_7_02	Standard	800	96.5%	<input type="checkbox"/>
DoubleBlank1	Double Blank		#N/A	<input type="checkbox"/>
DoubleBlank2	Double Blank		#N/A	<input type="checkbox"/>
DoubleBlank3	Double Blank		#N/A	<input type="checkbox"/>
QC_High_01	Quality Control	589	104.2%	<input type="checkbox"/>
QC_High_02	Quality Control	589	104.4%	<input type="checkbox"/>
QC_High_03	Quality Control	589	106%	<input type="checkbox"/>
QC_Low_01	Quality Control	121	116.6%	<input type="checkbox"/>
QC_Low_02	Quality Control	121	111.2%	<input type="checkbox"/>
QC_Low_03	Quality Control	121	111.1%	<input type="checkbox"/>
QC_Mid_01	Quality Control	346	121.4%	<input type="checkbox"/>
QC_Mid_02	Quality Control	346	125%	<input type="checkbox"/>
QC_Mid_03	Quality Control	346	121.6%	<input type="checkbox"/>
SPQC_01	Unknown		#N/A	<input type="checkbox"/>

The FDA Guidance on which this assay was based states that calibration points should have bias < 15% (accuracy between 85% and 115%) between the known concentration and the backcalculated concentration from the calibration curve, and from the Accuracy column we can see that Cal_5 does not meet that test. These replicates can be removed from consideration using the checkbox in the **Exclude from Calibration** column in the Document Grid, or by right-clicking on the outlier in the Calibration Curve view and selecting **Exclude from Calibration**.

- In the Document Grid, click on the checkbox in the **Exclude from Calibration** column for the **Cal5_01** replicate, then press the **down-arrow key**
- Repeat for **Cal5_02**

The Calibration Curve should now look like the image below. Note that the R-squared value improves from 0.97 to >0.99 by excluding the outlier Cal_5.



Now we can import the rest of our unknowns:

- On the **File** menu, select **Import** and click on **Results**.
- In the **Import Results** form, choose to import single-injection replicates. For best performance, be sure to select "Many" in the "Files to import simultaneously" control at the bottom of the form.
- Now click the **OK** button.

- You should see a file selection dialog with all the raw data folders in the tutorial folder. Select the as-yet-unloaded unknowns and click **OK**.

A convenient way to view quantitation data is to once again use the Document Grid, this time with the **Peptide Ratio Results** view.

- In the **Document Grid** window, click on the **View** control and select **Peptide Ratio Results**.
- Click on the **Replicate** column header and select **Sort Ascending**

The Document Grid should look something like this:


Document Grid: Peptide Ratio Results							
Peptide	Protein	Replicate	Peptide Peak Found Ratio	Peptide Retention Time	Ratio To Standard	Quantification	
Drug	DrugX	80_0_1_1_00_1...	1	2.68	0.6469	62.6089 uM	
Drug	DrugX	81_0_1_1_00_1...	1	2.68	0.4071	39.5176 uM	
Drug	DrugX	82_0_1_1_00_1...	1	2.68	0.5988	57.9825 uM	
Drug	DrugX	83_0_1_1_00_1...	1	2.68	0.3669	35.6509 uM	
Drug	DrugX	84_0_1_1_00_1...	1	2.68	0.3592	34.9036 uM	
Drug	DrugX	87_0_1_1_00_1...	1	2.68	0.2942	28.6499 uM	
Drug	DrugX	88_0_1_1_00_1...	1	2.68	0.531	51.4477 uM	
Drug	DrugX	89_0_1_1_00_1...	1	2.68	0.8523	82.3866 uM	
Drug	DrugX	90_0_1_1_00_1...	1	2.68	0.7145	69.1246 uM	
Drug	DrugX	91_0_1_1_00_1...	1	2.68	0.5668	54.897 uM	
Drug	DrugX	92_0_1_1_00_1...	1	2.68	0.7113	68.8121 uM	
Drug	DrugX	93_0_1_1_00_1...	1	2.68	0.5287	51.2276 uM	
Drug	DrugX	94_0_1_1_00_1...	1	2.68	0.773	74.7568 uM	
Drug	DrugX	95_0_1_1_00_1...	1	2.68	0.3779	36.7044 uM	
Drug	DrugX	96_0_1_1_00_1...	1	2.68	0.4428	42.957 uM	
Drug	DrugX	Blank_01	0.5	2.68	0	0.3157 uM	
Drug	DrugX	Blank_02	1	2.67	0.002	0.5048 uM	
Drug	DrugX	Blank_03	1	2.68	0.0011	0.4218 uM	
Drug	DrugX	Cal_1_01	1	2.68	0.105	10.4297 uM	
Drug	DrugX	Cal_1_02	1	2.68	0.0996	9.9036 uM	
Drug	DrugX	Cal_2_01	1	2.68	0.1913	18.7352 uM	
Drug	DrugX	Cal_2_02	1	2.68	0.2007	19.6453 uM	
Drug	DrugX	Cal_3_01	1	2.68	1.1231	108.4626 uM	
Drug	DrugX	Cal_3_02	1	2.68	1.0426	100.7119 uM	
Drug	DrugX	Cal_4_01	1	2.68	2.0557	198.2779 uM	
Drug	DrugX	Cal_4_02	1	2.68	1.9907	192.0188 uM	
Drug	DrugX	Cal_5_01	1	2.68	2.8707	276.759 uM	
Drug	DrugX	Cal_5_02	1	2.68	2.9928	288.5146 uM	
Drug	DrugX	Cal_6_01	1	2.68	6.3822	614.9043 uM	
Drug	DrugX	Cal_6_02	1	2.68	6.6078	636.6334 uM	
Drug	DrugX	Cal_7_01	1	2.68	8.3953	808.7618 uM	
Drug	DrugX	Cal_7_02	1	2.68	7.5365	726.0585 uM	
Drug	DrugX	DoubleBlank1	0	#N/A	#N/A	#N/A	
Drug	DrugX	DoubleBlank2	0	#N/A	#N/A	#N/A	
Drug	DrugX	DoubleBlank3	0	#N/A	#N/A	#N/A	
Drug	DrugX	QC_High_01	1	2.68	5.9928	577.41 uM	
Drug	DrugX	QC_High_02	1	2.68	6.0076	578.8339 uM	
Drug	DrugX	QC_High_03	1	2.68	6.0969	587.4288 uM	
Drug	DrugX	QC_Low_01	1	2.68	1.3809	133.2886 uM	
Drug	DrugX	QC_Low_02	1	2.68	1.3179	127.2255 uM	
Drug	DrugX	QC_Low_03	1	2.68	1.3159	127.0333 uM	
Drug	DrugX	QC_Mid_01	1	2.68	4.1029	395.4108 uM	
Drug	DrugX	QC_Mid_02	1	2.68	4.2251	407.1797 uM	
Drug	DrugX	QC_Mid_03	1	2.68	4.1094	396.0431 uM	

After removal of the two Cal_5 datapoints, further exploration of the data reveals that one of the Cal_7 points has accuracy <85% and should likely be removed. This will have little effect on the measurement of the samples, since there are no samples above the level of Cal_6, and only four samples have levels between Cal 4 and Cal 6.

To enable easier visualization of the dynamic range of the samples along the calibration curve:

- Right-Click in the calibration curve window and select **Log Plot**.

The calibration curve should look like this:



This presents a view which allows the user to easily observe that the samples mostly fall between Cal_2 (20 uM) and Cal_3 (100 uM), and well within the linear calibration range of the assay. The Quality Control samples (known unknowns, green diamonds above) all have accuracy measured between 85 and 115%, which meets the FDA Guidance criteria.

From here, the next steps would be to export the data for external statistical processing, or establish biological grouping within this document and utilize some of the statistical analysis tools or plugins within skyline. These options are covered in other tutorials.

Conclusion

In this tutorial, you have learned how to create a Skyline document that targets small molecules specified as precursor ion chemical formulas and adducts, and product ion m/z values. You imported a multi-replicate data set collected using LC-MS/MS on a triple quadrupole, and saw how many existing Skyline features created initially for targeted proteomics use can now be applied to small molecule data. Non-proteomic molecule support is still a relatively new feature area for Skyline. As such, you can expect it to continue improving rapidly.

Bibliography