


- Monitor instrument performance over time
- Maintain history and capture context
- For labs doing SRM, PRM, and DIA
 - Spike in controls
 - Run separate targeted QC samples (SRM or PRM)
- For labs doing DDA
 - Run separate targeted QC samples (PRM)
- Works with any instrument vendor/type
- Inspired by Statistical Process Control in Proteomics (SProCoP), a Skyline external tool

- Fully automated workflow with AutoQC
- New folder type within Panorama
- Expected usage pattern is one folder per instrument
- User interface focused on showing QC data
 - Levey-Jennings and Pareto plots
 - Not oriented around displaying specific Skyline documents

- Create Skyline document with the precursors/peptides/transitions to monitor
 - Typically a small number of peptides (10ish)
- Run QC sample on instrument
- Raw data file auto-loaded into Skyline document
- Data automatically pushed to Panorama
- View current data, identify issues
- Repeat

Levey-Jennings Plots

- Shows metric values over time
 - Retention time, peak area, etc
- One plot/series per peptide
- Date grouping and zooming




- Establish normal ranges for the assay
- Training data
 - Find period of stable performance
 - Defined by start and end dates
 - Calculate mean and standard deviation
- Future data is compared with training data
- Can have multiple guide sets per folder

- Capture additional context about the data
- Customizable, but typical annotations include:
 - Instrumentation changes
 - Reagent changes
 - Technician changes
- Overlaid on Levey-Jennings plots
- Helpful for correlating changes with causes

Pareto Plots

- Highlight most common QC issues
- One bar per metric
- New in upcoming 15.3 release



Demo

- Filter plots by other criteria
 - Lab tech
 - Study type
- Cross-instrument dashboard
 - Monitor multiple folders at once
 - Highlight likely problems
 - Ensure QC data is flowing into Panorama
- Improved statistical tools
- Automated notifications