






# Generating fit-for-purpose targeted PRM assays from DIA experiments

2023 Skyline User Meeting

Ariana E. Shannon, Yi Wang, Xin Gang, Amanda B. Hummon and Brian C. Searle


## CD8 Cells Kill Virally Infected Cells






# DIA validation of all targets for on-the-fly PRM


## 1. Chromatogram Library and Target Generation




## 2. Matrix-Matched curve for LOD/LOQ of all targets




## 3. Stability in short and long-term storage




## 4. Repeatability with and across digestions



# Dimethyl labelled reference peptides serve as a suitable background





RT deviation between 100% reference (unlabeled) and 10% reference injections




- Deviation of approximately 60 seconds with dimethyl labelled background

## Stability


### LOD/LOQ



### GRAB Stability Over 5 Freeze-Thaw Cycles




### GRAB Stability Over Time in 4C



### Peptide




- ASFRGDSGGPLVC[+57.0214635]K
- C[+57.0214635]IPHPDYNPK
- DQQPEAIC[+57.0214635]GGFLIR
- DREC[+57.0214635]ESYFK
- KVAAGIVSYGYK
- NVNVKPGDVC[+57.0214635]YVAGWGR
- RNVNVKPGDVC[+57.0214635]YVAGWGR
- TFSNDIMLLK
- TFSNDIMLLKLLK
- TNQIC[+57.0214635]AGDPK
- TNQIC[+57.0214635]AGDPKTK
- TQQVIPMVK
- VAAGIVSYGYK
- VAAGIVSYGYKDGSPPR
- VSSFLSWIK
- VSSFLSWIKK
- YNKTNQIC[+57.0214635]AGDPK
- YSNTLQVELTVQK
- YSNTLQVELTVQKDR

### Reproducibility (5 digestions x 5 days)



Figures of Merit quantified for 47,717 peptides identified in CD8<sup>+</sup> T cells after acutely stimulation. 13,398 peptides passed all validation assays.

# PRM with validated targets




25 proteins were detected at 0.3 ng level on the Exploris 480, equivalent to 10 T cells or 1 HeLa cell.


# Conclusions

- Dimethyl peptides can serve as a suitable background for matrix-matched calibration curves.
- Validating thousands of peptides from a DIA experiment allows for the generation of quick, preliminary PRM assays prior to acquiring stably isotopic labelled peptides.
- Mass spectrometry assays can be sensitive enough to compete with flow cytometry for monitoring protein abundance.

## Acknowledgements



Brian Searle  
Yi Wang  
Damien Wilburn



Amanda Hummon  
Emily Sekera  
Brian Fries

OSU  
Xin Gang  
Zihai Li  
Nojoon Soon